



NUTEC GLOBAL




General Surgical Instruments
Suction instruments


RHOTON SUCTION TIPS CONIC

Registration No.
878.4200


- The Rhoton Suction Tips Conic has excellent control of suction power, with the control hole on the anatomical face of the product.
- Thin micro-suction instrument.
- Atraumatic active tip.
- Wider view of the surgical site, even in the most restricted approaches.
- Single use.



Flexibility with **molding** and manual 360° rotational adjustment for multifunctionality in any surgical access:




- Easy connection with suction tube on stepped fitting:




- Tapered suction tube tip prevents clogging:








- Suction control in the controller hole:



- Comes with a disposable cleaning wire:



Dimensions

Description	Ø Tip		Working Length		Total Length		Color	Illustration
	MM	FR	MM	INCH	MM	INCH		
Rhoton Suction Tips Conic FR 3	1,0	3	80	3 1/8"	142	5 1/8"	Lilac	
			100	4"	162	6"		
			120	4 3/4"	182	6 1/2"		
			140	5 1/2"	202	7 1/4"		
			160	6 1/2"	222	8 3/4"		
			180	7"		9"		
Rhoton Suction Tips Conic FR 4	1,3	4	80	3 1/8"	142	5 1/8"	Green	
			100	4"	162	6"		
			120	4 3/4"	182	6 1/2"		
			140	5 1/2"	202	7 1/4"		
			160	6 1/2"	222	8 3/4"		
			180	7"		9"		
Rhoton Suction Tips Conic FR 5	1,7	5	80	3 1/8"	142	5 1/8"	Orange	
			100	4"	162	6"		
			120	4 3/4"	182	6 1/2"		
			140	5 1/2"	202	7 1/4"		
			160	6 1/2"	222	8 3/4"		
			180	7"		9"		
Rhoton Suction Tips Conic FR 6	2,0	6	80	3 1/8"	142	5 1/8"	Blue	
			100	4"	162	6"		
			120	4 3/4"	182	6 1/2"		
			140	5 1/2"	202	7 1/4"		
			160	6 1/2"	222	8 3/4"		
			180	7"		9"		
Rhoton Suction Tips Conic FR 7	2,3	7	80	3 1/8"	142	5 1/8"	Yellow	
			100	4"	162	6"		
			120	4 3/4"	182	6 1/2"		
			140	5 1/2"	202	7 1/4"		
			160	6 1/2"	222	8 3/4"		
			180	7"		9"		

The manufacturer recommends single use and not resterilizing the product.

Made in Brazil



NUTEC GLOBAL is committed to providing high-quality hospital medical products, incorporating technology and innovation to improve health and treatment effectiveness. Our commitment is to establish a prominent position in the surgical instruments and medical equipment market, driving solid and continuous growth. In addition, we are fully aligned with the requirements of RDC 665/2022, ensuring that our products meet the most rigorous quality and regulatory compliance standards.

Our mission translates into an unwavering dedication to patient well-being and healthcare professional satisfaction.

Through successful work and constant pursuit of improvement, we are committed to promoting significant advances in the area of medical products, always maintaining quality and innovation as our fundamental pillars. At NUTEC GLOBAL, we are determined to make a difference in people's health and lives.





NUTEC GLOBAL

NUTEC GLOBAL INDÚSTRIA E COMÉRCIO LTDA



Av. 52 A, 771 , Vila São José,
Rio Claro, São Paulo, Brasil.

comercial@nutecglobal.tech

+ 55 019 3536-3535

www.nutecglobal.com.br

CAT.0000.003 | Revisão 00