




FRAZIER SUCTION TIP CONIC

FRAZIER SUCTION


FRAZIER SUCTION
TIP CONIC



Cleaning thread

Registration No.
878.4200


Technical Specifications

Description	Ø Tip		Working Length		Total Length		Color	Illustration
	MM	FR	MM	INCH	MM	INCH		
Frazier Suction Tip Conic	3,3	10	80	3 1/8"	142	5 1/8"	Pink	
			100	4"	162	6"		
			120	4 3/4"	182	6 1/2"		
			140	5 1/2"	202	7 1/4"		
			160	6 1/2"	222	8 3/4"		
			180	7"		9"		

Flexibility with **molding** and manual 360° rotational adjustment for multifunctionality in any type of surgery




Moldable tubes



- Sterile Product - Sterilized by Ethylene Oxide.
- Shelf life of 2 years.
- The manufacturer recommends Single Use and Do Not Resterilize the product.


- Suction control in the controller hole




- Comes with a disposable cleaning cord



- Easy connection with suction tube on stepped fitting



- Tapered tube tip prevents clogging



- Wider view of the operating site, even in the most restricted approaches.
- Single-use
- Atruamatic Active Tip
- Moldable instrument

Nutec Global Indústria e Comércio LTDA

Avenida 52A, nº 771 - Vila Alemã - Rio Claro/ SP - CNPJ: 07.575.813/0001-59 |
+55 (19) 3536-3535 | nutecglobal.com.br | comercial@nutecglobal.tech



CAT.0000.004 | Revisão 00